## GATE INSTALLATION FRAMING GATE SILL FLANGE LEFT HAND GATE INSTALLATION SHOWN SHOULD BE PLACED FLUSH WITH UNDERNEATH FRAMING-NOTES: GATE OPENS 1. ALL FRAMING MUST BE CAPABLE OF WITHSTANDING A 200 LBS CONCENTRATED LOAD IN ALL DIRECTIONS. THIS 2. AN ALTERNATIVE METHOD TO ACHIEVE A FLUSH FINISH BETWEEN THE GATE SILL FLANGE AND FASCIA BOARD IS TO ROUTER A POCKET INTO THE THE AREA COVERED BY THE GATE SILL FLANGE. THIS POCKET SHOULD BE 3/16" DEEP. 3. REFER TO THE TECHNICAL DRAWINGS FOR ADDITIONAL INSTALLATION DIMENSIONS. 44 1/4" (36" WIDE GATE) 50 1/4" (42" WIDE GATE) DETAIL B SCALE 1:2 -4" X 4" 4" X 4"--3/16" 42 1/4" OPEN SPACE TO ALLOW GATE TO BE INSTALL MEASURED FROM 4 7/8" TOP OF FINISH FLOOR TO FASCIA BOARD -3/8"MINIMUM 2"x4"s TO BE PLACED BACK TO ALLOW GATE ANGLE AND FASCIA BOARD TO BE FLUSH WITH 4"X4" POSTS, CUT WOOD AS NEEDED -2" X 4" DOUBLE FRAMED 53", 75",123",147" TO FIT SPECIFIC APPLICATION. DETAIL A OR 171" SCALE 1:2 (MODEL SPECIFIC) FASCIA BOARD FASCIA BOARD CUTAWAY MAX RECOMMEND 3/8" MINIMUM THICK TO SHOW FRAMING FLOOR PANEL 2" X 4" TO BE CONCRETE SLAB **FLOOR** ANCHORED TO NOTE: CONCRETE SLAB-THE LENGTH & WIDTH OF THE SLAB MAY NEED TO BE INCREASED OVER THE MINIMUM SLAB DIMENSIONS SHOWN ON THE TECHNICAL DRAWINGS SECTION A-A

ILS-00804 SHEET 1 OF 1 REV. 4 (6825)(1/04/19)(BG)

