



Specifically designed to withstand harsh environments, Bruno's Outdoor Elite stairlift gives you the ability to enjoy the outside of your home safely. Whether you need help getting up porch stairs or down steps to the pier, the Outdoor Elite reliably takes you there with a smooth, safe ride.

- Weight capacity: 400 lb (181 kg)
- Performance tested for 0° F to 125° F (-18° C to 52° C)
- Durable outdoor paint protects all steel surfaces
- Anodized aluminum rail with covered gear rack
- Marine-grade vinyl padded seat and armrest
- Reliable performance even in a power outage
- Generous size seat with multiple height adjustments
- Onboard audio/visual diagnostics for easy service
- Tight rail installation: installs to within 7 in (178 mm) of the wall
- Extended rail lengths available
- Gold warranty: five year major components; two year parts

See reverse side for additional features.

ELITE OUTDOOR STRAIGHT FEATURES



TECHNICAL SPECIFICATIONS

Seat Pad Size: 17.5 in wide x 15 in deep
(445 mm wide x 381 mm deep)

Width Between Arms: Adjustable to 21.75 in
(552 mm) between arms

Footrest Size: 10 in long x 9 in wide
(254 mm long x 229 mm wide)

Learn more at www.bruno.com



Footrest/carriage safety sensors stop the unit when obstruction encountered.



Offset swivel seat makes getting on/off easy.



Flip up arms, seat and footrest create plenty of extra space on steps.



Weather-resistant cover provides durable protection.



Two wireless call/send controls.



Bruno Independent Living Aids, Inc.
1780 Executive Drive, Oconomowoc, WI 53066
www.bruno.com

© Bruno Independent Living Aids, Inc. 2017
2017-03-23-Elite-Outdoor-Straight-.indd • MKT-SL-4