



Easily access all levels of your home

High-quality American craftsmanship is the hallmark of custom-made Bruno Elite curved rail stairlifts. Designed to the exact specifications of your staircase, the Elite melds into every curve. The plush chair properly reflects the beauty of your home.

Let Bruno's Elite safely connect you to all levels of your home again.

- Made-in-the-USA quality
- Custom-designed rail built for your stairs
- Smooth, quiet ride

Learn more at www.bruno.com.

Even a Legend Can Use a Lift



“We completed a sports room in our home to enjoy the football season with our family just before Bart had his strokes. Afterwards, it was not possible to use this area because Bart was unable to negotiate the stairs. A Bruno stairlift allows us to enjoy our entire home again.”

Bart & Cherry Starr

“Green Bay Packers Hall of Fame Quarterback”




Luxury styling and comfort.


ELITE CURVE OPTIONS



Top or bottom park position extends rail away from the staircase.



Power folding footrest automatically flips up/down when seat is raised/lowered.



Power swivel seat for effortless exit. Controlled on chair arm or wireless call/send.



Mid-park and charge station for staircases with middle landings.



Offset swivel seat makes it easy to get on/off. Add the power seat option to make it effortless.

ELITE CURVE FEATURES

- Weight capacity: 400 lb (181 kg)
- Plush, generous-size seat with multiple height adjustments
- Adjustable width armrests to 21.75 in (552 mm) between arms
- Offset swivel seat makes getting on/off safe and easy
- Armrest rocker switch for easy operation
- Painted rail with covered gear rack gives a sleek appearance
- Flip up arms, seat, and footrest create plenty of extra space on steps
- Two wireless call/sends allow unit to be controlled remotely
- Two 12v batteries ensure performance even in a power outage
- Soft start/stop for a smooth ride from start to finish
- Adjustable footrest height for maximum comfort
- Onboard audio/visual diagnostics for easy service
- Obstruction sensors and seat belt ensure safety
- Gold warranty: five-year major components; two-year parts



Custom crafted for a precision fit around curves.

UPHOLSTERY

Choose the Elite's standard tan vinyl or one of many optional colors and fabrics. Ask your local dealer for samples to see what best matches your home.





TECHNICAL SPECIFICATIONS

Standard Seat Size: 17.5 in wide x 15 in deep
(445 mm wide x 381 mm deep)

Larger Seat Option: 21.5 in wide x 16 in deep
(546 mm wide x 406 mm deep)

Standard Footrest Size: 10 in long x 9 in wide
(254 mm long x 229 mm wide)

Larger Footrest Option: 12 in long x 11 in wide
(305 mm long x 279 mm wide)

About Bruno

Bruno Independent Living Aids has helped improve the lives of people with limited mobility for over 30 years. Engineered and manufactured in the USA, Bruno stairlifts and vertical platform lifts help people better access homes, public buildings and businesses, and Bruno's full line of scooter/powerchair lifts help people easily transport their mobility devices. In addition, Bruno's Valet Signature Seating is the world leader in revolutionary turning automotive seats to assist people in getting in and out of their vehicles more easily. A family-owned business, Bruno takes pride in their unwavering focus on quality and safety and their position as industry innovation leader.

Proudly made in the USA.



Bruno Independent Living Aids, Inc.
1780 Executive Drive, Oconomowoc, WI 53066

www.bruno.com