



## TECHNICAL SPECIFICATIONS

**Lifting Capacity:** 100 lb (45 kg)

**Vehicle Type:** CUVs, hatchbacks, minivans, pickups, sedans, SUVs, full-size vans

**Hitch Type:** I, II or III  
(receiver hitch required)

Learn more at [www.bruno.com](http://www.bruno.com)

Lift and stow a folding manual wheelchair with Bruno's Back-Saver. Roll the wheelchair on the platform, touch the button, and your wheelchair will be raised into place and automatically secured for travel. Designed to handle 16" and 22" wide folding wheelchairs with 22" or larger wheels.

- Power up/down platform
- Self-locking securement
- Manual backup for peace of mind
- Potential to reinstall in another applicable vehicle
- Eligible for auto manufacturer rebate\*
- 3-year limited warranty

\* Contact Bruno dealer for details.

## BACK-SAVER



## BACK-SAVER OPTIONS



- Vinyl cover (not shown)
- License plate holder (not shown)
- Swing-Away for easy access to hatch/cargo space (left)



Bruno Independent Living Aids, Inc.  
1780 Executive Drive, Oconomowoc, WI 53066  
[www.bruno.com](http://www.bruno.com)

© Bruno Independent Living Aids, Inc. 2019  
2019-04-10 Back-Saver.indd • MKT-L-9