



✓ Verified order

"The stair lift is great, very solid and looks good too."

L. Uhring



TECHNICAL SPECIFICATIONS

Standard Seat Size: 17.5 in wide x 15 in deep
(445 mm wide x 381 mm deep)

Larger Seat Option: 21.5 in wide x 16 in deep
(546 mm wide x 406 mm deep)

Standard Footrest Size: 10 in long x 9 in wide
(254 mm long x 229 mm wide)

Larger Footrest Option: 12 in long x 11 in wide
(305 mm long x 279 mm wide)

Vertical Rail:

Installs within 5 in (127 mm) of the wall

INDUSTRY-LEADING WARRANTY

Lifetime Warranty on Major Components -- motor, gearbox, drivetrain and rail -- for the lifetime of the Bruno Stairlift as long as the original purchaser owns the stairlift.

2-Year Warranty on all components after date of installation by an Authorized Bruno Dealer.

30-Days Labor Coverage after date of installation by an Authorized Bruno Dealer for any manufacturing defect.

About Bruno

Family-owned Bruno Independent Living Aids has improved the lives of people with limited mobility for over 35 years. Engineered and manufactured in the USA, Bruno makes stairlifts and vertical platform lifts to help people better access homes, public buildings and businesses and scooter/powerchair lifts to assist people in easily transporting their mobility devices.

Proudly made in the USA.



Bruno Independent Living Aids, Inc.
1780 Executive Drive, Oconomowoc, WI 53066
www.bruno.com

© Bruno Independent Living Aids, Inc. 2020
2020-03-16-Elite-6-PageV8.indd • MKT-SL-2



CONSUMERAFFAIRS



BRING HOME **COMFORT & ACCESS**

www.bruno.com





COMBINE QUALITY, COMFORT AND STYLE

When you want the highest quality stairlift combined with a stylized appearance, you want the Bruno Elite. The Elite offers worry-free access to all levels of your home. And the Elite's stylized design blends with your home's environment.

Let Bruno's Elite safely connect you to all levels of your home again.

Made-in-the-USA
400 lb Lift Capacity (181 kg)



Verified order

"The fact they produce their products in their own American factory, and the favorable reviews impressed me. Seeing the quality of the chairlift in our home cemented my belief we made the perfect decision with Bruno."

Tim P.

HIGH-PERFORMANCE. EXCEPTIONAL RIDE.

- **Luxury Styling** - Plush, generous-size seat with multiple height adjustments, adjustable arm widths and footrest height for maximum comfort.
- **Space Saving** - Narrow vertical rail preserves maximum stairway space. Fold arms, seat and footrest to save even more space.
- **Simple Operation** - Control on armrest activates stairlift. Two wireless remote controls also included.
- **Easy On/Off** - Seat swivels up to 90 degrees at top landing, extends and locks in place for safe exit away from steps.
- **Sleek & Clean** - Covered gear rack on vertical rail gives clean appearance and reduces exposure to dirt, grease.
- **Safe** - Retractable seat belt and obstruction sensors ensure safety.
- **Dependable** - Battery-powered performance even in a power outage. Stairlift continually charging along entire rail.





Offset swivel seat makes it easy to get on/off.

★★★★★ Verified order

“Our Bruno chair was a wise and wonderful investment into our home to ensure that my husband and I can remain in our home for years to come.”

J. Therrien



Compact design for open space on stairs.

ELITE OPTIONS


- Multiple colors and fabrics



Power folding footrest automatically flips up/down when seat is raised/ lowered.



Power swivel seat for effortless exit. Controlled on chair arm or wireless call/send.



Power or manual folding rails for narrow hallway or when doorway is at the bottom of stairs.



Larger seat and/or footrest for increased comfort.